How to solve problems? Troubleshooting Guide

Problem	Discomfort	Solution
Air leaks with nasal mask	Thirst and mouth dryness	Check for mouth leaks
(1st choice)		Use chin strap
		Use facial mask
Humidification	Nasal dryness, congestion, and	Add heated humidifier
	rhinitis, nasal bleeding	Increase humidification
		Add heated circuit
Sleep	Sleep fragmentation	Check for air leaks or upper airway obstruction
Asynchronies	Aerophagia, eructation, flatulence,	Reconsider inspiratory trigger synchrony and backup RR
	and abdominal discomfort	Decrease IPAP
Deventilation dyspnea	Morning breathlessness	Decrease IPAP
Inspiratory trigger	Difficulty to trigger the breath	Increase inspiratory trigger sensitivity
asynchrony		Increase EPAP if intrinsic PEEP suspected
Inspiratory trigger	Auto trigger	Check for non-intentional leaks
asynchrony		Decrease inspiratory trigger sensitivity
Flow asynchrony	Air arrives too quickly	Prolong rise time
		Decrease IPAP
Flow asynchrony	Not enough air	Decrease rise time
		Increase IPAP
Cycling asynchrony	Breath too short	Prolong expiratory trigger and maximum inspiratory time
Cycling asynchrony	Double triggering	Prolong expiratory trigger and maximum inspiratory time
		Increase IPAP
Cycling asynchrony	Breath too long	Shorten expiratory trigger and minimum inspiratory time
Rate asynchrony	No time to exhale	Decrease backup RR, Check Insp Time?
		EDUCATION by Breas